

# Surgicare<sup>®</sup>

## syntho premium

Latexfree and Powderfree Nitrile Medical Examination Gloves

# Surgicare<sup>®</sup>

## syntho

Latexfree and Powderfree Nitrile Medical Examination Gloves

Specially designed Latex Free and Powder free Nitrile Examination gloves, to suit every Healthcare Professional. It is a glove that provides excellent barrier protection for the healthcare worker against potentially harmful bodily fluids and hazardous drugs.

### SPECIAL FEATURES

- ▶ Latex free and Powder free
- ▶ Ambidextrous for quick and easy donning
- ▶ Finger textured for additional grip
- ▶ Beaded cuff for easy donning
- ▶ Length from 230 mm
- ▶ Weight of 3.5 gms. (3.2 mil) and 5.5 gms. (5 mil)
- ▶ ASTM : D-6319-2019.
- ▶ Certified to meet PPE standards EN 374, EN 420, Cat III
- ▶ Device classification - Class 1

### Nitrile Exam gloves is the perfect choice for

- ▶ Critical Care Units & ICU's.
- ▶ Dental Workers.
- ▶ Healthcare Workers administering potentially hazardous drugs and for clean up of spills and handling waste of treated patients.
- ▶ Oncology Units and Cancer Research Institutes handling cytotoxic drugs (5mil Nitrile Gloves).
- ▶ Ophthalmology units and Pathology units.
- ▶ For use in the food industry and hospitality sector.

### QUALITY STANDARDS

- ▶ Manufactured in accordance to ISO 9001: 2015 and ISO 13485: 2016 Quality System Standards.
- ▶ Bio-compatible as measured by primary skin sensitization test as per ISO 10993 -2010
- ▶ Tested for use with 9 Chemotherapy Drugs as per ASTM D 6978 - 05 (for 5 mil Nitrile Exam gloves - Surgicare Syntho Premium)
- ▶ Food grade certified.

### QUALITY CONTROL

- ▶ 100% Inspected
- ▶ AQL - 1.5

[www.kanamlatex.com](http://www.kanamlatex.com)



# Kanam

## KANAM LATEX INDUSTRIES PVT. LTD.

The trusted name the world over

Ooppotttil Buildings, K.K. Road  
Kottayam-686 001, Kerala, India  
Tel: 0481-2560108, 2563513  
Fax: 0481-2563614

E-mail: [customercare@kanamlatex.com](mailto:customercare@kanamlatex.com)  
[ho@kanamlatex.com](mailto:ho@kanamlatex.com)

## Product Information

### TYPE

Powder Free Nitrile Examination Gloves.

### PRIMARY MATERIAL

Nitrile

### DESIGN

Ambidextrous, Beaded cuff and textured.

### POWDER FREE RESIDUE

<2mg/glove

**Tested for use with Chemotherapy Drugs: Per ASTM D 6978-05**  
(Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs).

**The following chemotherapy drugs and concentration had no breakthrough detected upto 240 minutes**

(for 5 mil Nitrile Exam gloves - Surgicare Syntho Premium):

- ▶ **Cisplatin**, 1.0 mg/ml (1,000ppm)
- ▶ **Cyclophosphamide (Cytosan)**, 20 mg/ml (20,000ppm)
- ▶ **Dacarbazine (DTIC)**, 10.0 mg/ml (10,000ppm)
- ▶ **Doxorubicin Hydrochloride**, 2.0 mg/ml (2,000ppm)
- ▶ **Etoposide (Toposar)**, 20.0 mg/ml (20,000ppm)
- ▶ **Fluorouracil**, 50.0 mg/ml (50,000ppm)
- ▶ **Paclitaxel (Taxol)**, 6.0 mg/ml (6,000ppm)

**The following chemotherapy drugs and concentration showed breakthrough detected in less than 30 minutes** (for 5 mil Nitrile Exam gloves - Surgicare Syntho Premium):

- ▶ **Carmustine** (3.3 mg/ml): 1.8 minutes
- ▶ **ThioTEPA** (10.0 mg/ml): 1.6 minutes.

**Warning:** Not for use with: Carmustine, ThioTEPA

## PHYSICAL DIMENSIONS

Product	Length	Palm Width (mm)			
		S	M	L	XL
Surgicare Syntho (3.2 mil)	230	84 - 85	96 - 97	106 - 107	112 - 113
Surgicare Syntho Premium (5 mil)	230	84 - 85	96 - 97	106 - 107	112 - 113

## PHYSICAL PROPERTIES

Characteristics	Before Ageing	After Ageing
Tensile Strength (Mpa) min.	14	14
Ultimate Elongation (%) min.	500	400

## PACKAGING

Material	Syntho Premium	Syntho
Duplex Box (Quantity)	100 pcs.	200 pcs.
Master Carton (Quantity)	1000 pcs.	2000 pcs.

KIPUB/21/Rev# 02



**Surgicare<sup>®</sup>**  
**syntho premium**  
Latexfree and Powderfree Nitrile Medical Examination Gloves

**Surgicare<sup>®</sup>**  
**syntho**  
Latexfree and Powderfree Nitrile Medical Examination Gloves